

Product Datasheet

Human recombinant ACE2-Fc (orb1108525)

Description	Human recombinant ACE2-Fc
Preservatives	PBS pH 7.4
Storage	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles. Shipping on dry ice.
Note	For research use only
Source	CHO-based cell line (expressed by QMCF Technology)
MW	112,14 kDa; Protein contains amino acids 18-740 and human IgG1 Fc in C-terminus.
Uniprot ID	Q9BYF1
Expiration Date	6 months from date of receipt.