

## Product Datasheet

### SARS-CoV-2 Spike S2 (orb1108527)

<b>Description</b>	SARS-CoV-2 Spike S2
<b>Preservatives</b>	PBS pH 7.4
<b>Storage</b>	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles. Shipping on dry ice.
<b>Note</b>	For research use only
<b>Source</b>	CHO-based cell line (expressed by QMCF Technology)
<b>MW</b>	58,82 kDa; Protein contains amino acids 693-1218 and His-6 tag at C-terminus.
<b>Uniprot ID</b>	<a href="#">P0DTC2</a>
<b>Expiration Date</b>	6 months from date of receipt.