

Product Datasheet

SARS-CoV-2 S1 RBD (orb1108529)

Description	SARS-CoV-2 S1 RBD
Preservatives	PBS pH 7.4
Storage	Store at -70°C upon receipt. Recommended to aliquot into smaller quantities. Avoid repeated freeze-thaw cycles. Shipping on dry ice.
Note	For research use only
Source	CHO-based cell line (expressed by QMCF Technology)
MW	23 kDa; Protein contains amino acids 331-524, two extra amino acids (YS) in N-terminus and His-6 tag at C-terminus.
Uniprot ID	P0DTC2
Expiration Date	6 months from date of receipt.