

Product Datasheet

Monoclonal antibody to Procalcitonin, clone 9D8 (orb1108561)

Description	Monoclonal antibody to Procalcitonin, clone 9D8
Species/Host	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant hPCT protein produced by CHO-based cell line (expressed by QMCF Technology).
Preservatives	PBS pH 7.4
Storage	The product is stable at 4 °C for several days. Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	human IgG1-kappa
Clonality	Recombinant
Uniprot ID	P01258
Expiration Date	12 months from date of receipt.