

Product Datasheet

Monoclonal antibody to SARS-CoV-2 Nucleocapsid protein, clone 11H11 (orb1108568)

Description	Monoclonal antibody to SARS-CoV-2 Nucleocapsid protein, clone 11H11
Species/Host	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	SARS-CoV-2 virus (PBMCs isolated from patients recovered from SARS-CoV-2 infection, antibodies discovered by using HybriFree Technology).
Preservatives	PBS pH 7.4
Storage	The product is stable at 4 °C for several days. Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	human IgG1-lambda
Clonality	Recombinant
Uniprot ID	NA
Expiration Date	12 months from date of receipt.