

Product Datasheet

Monoclonal antibody to SARS-CoV-2 Spike RBD, clone 79C3 (orb1108572)

Description	Monoclonal antibody to SARS-CoV-2 Spike RBD, clone 79C3
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant SARS-CoV-2 RBD protein
Preservatives	PBS pH 7.4
Storage	The product is stable at 4 °C for several days. Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	rabbit IgG-kappa
Clonality	Recombinant
Uniprot ID	P0DTC2
Dilution Range	Indirect ELISA 0.01 – 0.02 ng/ml
Expiration Date	12 months from date of receipt.