

Product Datasheet

Monoclonal antibody to SARS-CoV-2 Nucleocapsid protein, clone 82C3 (orb1108579)

Description	Monoclonal antibody to SARS-CoV-2 Nucleocapsid protein, clone 82C3
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid protein
Preservatives	PBS pH 7.4
Storage	The product is stable at 4 °C for several days. Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	rabbit IgG-lambda
Clonality	Recombinant
Uniprot ID	P0DTC9
Dilution Range	Indirect ELISA 0.02 – 0.04 ng/ml
Expiration Date	12 months from date of receipt.