

Product Datasheet

Monoclonal antibody to SARS-CoV-2 Spike RBD, clone 36A7 (orb1108581)

Description	Monoclonal antibody to SARS-CoV-2 Spike RBD, clone 36A7
Species/Host	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 RBD2 protein
Preservatives	PBS pH 7.4
Storage	The product is stable at 4 °C for several days. Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	human IgG1-lambda
Clonality	Recombinant
Uniprot ID	P0DTC2
Dilution Range	Indirect ELISA 0.04 – 0.08 ng/ml
Expiration Date	12 months from date of receipt.