

Product Datasheet

Monoclonal antibody to SARS-CoV-2 Nucleocapsid protein, clone 16C11 (orb1108592)

Description	Monoclonal antibody to SARS-CoV-2 Nucleocapsid protein, clone 16C11
Species/Host	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid protein
Preservatives	PBS pH 7.4
Storage	The product is stable at 4 °C for several days. Store at -20 °C to -70 °C long term. Avoid multiple freeze-thaw cycles.
Note	For research use only
Isotype	human IgG1-lambda
Clonality	Recombinant
Uniprot ID	P0DTC9
Dilution Range	ELISA: 0.08 - 0.16 µg/ml WB: 0.5 - 1 µg/ml
Expiration Date	12 months from date of receipt.