

Product Datasheet

Chicken Anti-SARS-CoV-2 NP polyclonal antibody (orb1108593)

Description	Chicken Anti-SARS-CoV-2 NP polyclonal antibody
Species/Host	Gallus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 Nucleocapsid protein
Preservatives	Concentrated ammonium sulphate in PBS pH 7.4
Storage	Store at +4°C upon receipt. As product is (NH ₄) ₂ SO ₄ precipitate, mix well by pipetting or vortexing prior use.
Note	For research use only
Isotype	IgY
Clonality	Polyclonal
Uniprot ID	P0DTC9
Dilution Range	ELISA: 0,02-0,1 µg/ml WB: 10 µg/ml
Expiration Date	12 months from date of receipt.