

Product Datasheet

Mouse monoclonal antibody to HIV-1 Nef (orb1108671)

Description	Mouse monoclonal antibody to HIV-1 Nef
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	Recombinant Nef protein of HIV-1 BRU (purified from E. coli)
Preservatives	PBS pH 7.4, with 0.1% sodium azide
Storage	Store at -20...-70°C upon receipt. Divide antibody into aliquots prior usage. Avoid multiple freeze-thaw cycles
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Uniprot ID
Dilution Range	ELISA 0.027 to 0.8 µg/ml; WB 0.8 µg/ml
Expiration Date	12 months from date of receipt.