

Product Datasheet

Crry antibody (orb1108800)

Description	Mouse monoclonal antibody to Crry
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC, Functional Studies, IA, IF, IHC-Fr, WB
Form/Appearance	1 ml (100 µg/ml) 0.2 µm filtered protein G purified antibody solution in PBS, containing 0.1% bovine serum albumin.
Storage	Product should be stored at 4°C. Under recommended storage conditions, product is stable for at least one year.
Note	For research use only
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.