

Product Datasheet

CD83 antibody (orb1108827)

Description	Mouse monoclonal antibody to CD83
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IA, WB
Form/Appearance	1 ml (100 µg/ml) 0.2 um filtered antibody solution in PBS, containing 0.1% bovine serum albumin and 0.02% sodium azide.
Storage	Product should be stored at 4°C. Under recommended storage conditions, product is stable for at least one year.
Note	For research use only
Application notes	FC: Antibody ELBE-1 stains the extracellular domain of CD83. (Ref.1) W: A non-reduced and/or reduced sample treatment and SDS-Page was used. The band size is 40kDa kDa (Ref.1).
Isotype	Mouse IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.