

Product Datasheet

TLR2 antibody (Biotin) (orb1108996)

Description	Mouse monoclonal antibody against TLR2 conjugated to Biotin
Species/Host	Mouse
Reactivity	Human
Conjugation	Biotin
Tested Applications	FC, IA, WB
Form/Appearance	0.5 ml (100 µg/ml) 0.2 µm filtered biotinylated antibody solution in PBS, containing 0.02% sodium azide and 0.1% bovine serum albumin.
Storage	Product should be stored at 4degC. Under recommended storage conditions, product is stable for one year.
Note	For research use only
Isotype	Mouse IgG2a
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.