

Product Datasheet

Polymyxin B antibody (orb1109015)

Description	Mouse monoclonal antibody to Polymyxin B
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IA, WB
Form/Appearance	1 ml > 200 µg/ml) culture medium with a low endotoxin level containing 0.02% sodium azide.
Storage	Product should be stored at -20degC. Under recommended storage conditions, product is stable for one year.
Note	For research use only
Isotype	Mouse IgM
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.