

Product Datasheet

Crry/p65 antibody (orb1109086)

Description	Mouse monoclonal antibody to Crry/p65
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IA, IF, WB
Form/Appearance	1 ml (100 µg/ml) 0.2 µm filtered antibody solution in PBS, containing 0.1% bovine serum albumin and 0.02% sodium azide.
Storage	Product should be stored at 4°C. Under recommended storage conditions, product is stable for one year.
Note	For research use only
Application notes	IA: HM1117 can be used as detection antibody. IF: Positive signal was obtained with non-stimulated and RANKL-stimulated RAW264.7 cells. Concentration HM1117 used was 4µg/ml. W: A reduced sample treatment was used. the expected band size is ~50 kDa.
Isotype	Rat IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.