

Product Datasheet

Nop5p/Nop6p antibody (orb1109255)

Description	Mouse monoclonal antibody to Nop5p/Nop6p
Species/Host	Mouse
Reactivity	Other
Conjugation	Unconjugated
Tested Applications	ICC
Immunogen	Nop5p (a.k.a. Nop58p) and Nop6p (a.k.a. Nop56p). Both are localized in the nucleolus.
Preservatives	PBS, pH 7.4 with 0.05% sodium azide.
Form/Appearance	Liquid
Concentration	1ug/ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	Immunofluorescence on yeast cells: 1:500 to 1:5,000.
Clonality	Monoclonal
Dilution Range	Immunofluorescence on yeast cells: 1:500 to 1:5,000.
Expiration Date	12 months from date of receipt.