



Product Datasheet

N Cadherin antibody (orb11100)


Description	N Cadherin antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human N Cadherin. Please contact us for the exact immunogen sequence. The peptide is available as orb12897.
Target	N Cadherin
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	100 kDa
Uniprot ID	P15116, Q9Z1Y3, P19022
NCBI	001295105, 001308176
Entrez	1000
Dilution Range	IF/ICC: 1:50-500, IHC-P: 1:50-500, WB: 1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of rat heart (lane...)



IHC-P image of rat brain tissue using an...



Immunohistochemical staining of pig kidn...