



Product Datasheet

NCX1 antibody (orb11105)


Description	NCX1 antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human NCX1 (801-900/971 aa)
Target	NCX1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	106 kDa
Uniprot ID	P32418
Dilution Range	WB=1:500-2000, Flow-Cyt=0.2ug/test, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



U-937 cells were fixed with 4% PFA for 1...



permeabilized with 20% PBST for 20 min a...



Immunohistochemical staining of paraffin...