



Product Datasheet

NCX3 antibody (orb11106)


Description	NCX3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human NCX3 (641-740/927 aa)
Target	NCX3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100 kDa
Uniprot ID	P57103
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of rat brain...



Immunohistochemical analysis of paraffin...



Western blot analysis using NCX3 antibod...