

Product Datasheet

NeuN antibody (orb11112)

Description

NeuN antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

ELISA, IHC-P, WB

Applications

Immunohistochemical
analysis of paraffin...

Immunogen

KLH conjugated synthetic peptide derived from human NeuN. Please contact us for the exact immunogen sequence. The peptide is available as orb374739.

Target

NeuN

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl):
0.5 mg/ml

IHC-P of mouse brain
tissue (NeuN antibo...

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

IHC-P image of
mouse brain tissue
using ...

MW

34 kDa

Uniprot ID

[Q8BIF2](#), [A6NFN3](#)

NCBI

[21](#), [001024931](#), [001020102](#)

Entrez

[146713](#)

Dilution Range

IHC-P: 1:100-500 (based on 0.5 mg/ml), WB: 1:100-400

Expiration Date

12 months from date of receipt.