



Product Datasheet

NGF beta antibody (orb11126)


Description	NGF beta antibody
Species/Host	Rabbit
Reactivity	Canine, Human, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human NGF beta (151-220/241 aa)
Target	NGFB
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	13/32 kDa
Uniprot ID	P01138
Dilution Range	IHC-P=1:100-500, ICC=1:100, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



IHC-P of mouse hippocampus tissue (orb11...



Western blot analysis of rat brain lysat...