



Product Datasheet

NIK antibody (orb11130)


Description	NIK antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human NIK (901-947/947 aa)
Target	MAP3K14/NFkB Inducing Kinase
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	104 kDa
Uniprot ID	Q99558
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, Flow-Cyt=3ug/Test, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of Mouse Kidney ly...



Immunohistochemical analysis of human co...



Flow cytometric analysis of A431 cells U...