

Product Datasheet

Na⁺/K⁺-ATPase alpha 2 antibody [Out of stock] (orb1113117)

Description	Rabbit polyclonal antibody to Na ⁺ /K ⁺ -ATPase alpha 2
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from the C-terminal region of human Na ⁺ /K ⁺ -ATPase α2.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500 - 1:2000; ELISA 1:20000;
Expiration Date	12 months from date of receipt.