



Product Datasheet

LGR5 antibody (orb11148)


Description	Rabbit polyclonal antibody to LGR5
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant fusion protein of human LGR5/GPR49(NP_003658.1).
Preservatives	PBS with 0.01% thiomersal, 50% glycerol, pH 7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	100kDa
Uniprot ID	075473
Dilution Range	WB: 1:1000 - 1:2000, IF/ICC: 1:50 - 1:200
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of PC12 cell...



Immunofluorescence analysis of NIH/3T3 C...



Western blot analysis of extracts of var...