



Product Datasheet

Cd304/Bdca-4/Nrp-1 antibody (orb11167)


Description	Cd304/Bdca-4/Nrp-1 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human neuropilin-1 (851-923/923 aa)
Target	Neuropilin 1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	103 kDa
Uniprot ID	O14786
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



IHC-P of neural cilia of human colon can...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of formalin...