

Product Datasheet

Synaptotagmin 1/2 (phospho-S309/306) antibody [Out of stock] (orb1116707)

Description	Rabbit polyclonal antibody to Synaptotagmin 1/2 (phospho-S309/306)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from human Synaptotagmin 1/2 around the phosphorylation site of S309/306.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500 - 1:2000; IHC 1:100 - 1:300; ELISA 1:10000;
Expiration Date	12 months from date of receipt.