

Product Datasheet

Tubulin alpha 1/3/4 (phospho-Y272) antibody [Out of stock] (orb1116791)

Description	Rabbit polyclonal antibody to Tubulin alpha 1/3/4 (phospho-Y272)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from human Tubulin α 1/3/4 around the phosphorylation site of Y272.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500 - 1:2000; ELISA 1:40000;
Expiration Date	12 months from date of receipt.