

## Product Datasheet

### Cleaved-Caspase-9 p35 (D315) antibody [Out of stock] (orb1117439)

|                            |                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to Cleaved-Caspase-9 p35 (D315)                                                            |
| <b>Species/Host</b>        | Rabbit                                                                                                                |
| <b>Reactivity</b>          | Human                                                                                                                 |
| <b>Conjugation</b>         | Unconjugated                                                                                                          |
| <b>Tested Applications</b> | ELISA, IF, IHC-P, WB                                                                                                  |
| <b>Immunogen</b>           | Synthesized peptide derived from the Internal region of human Caspase-9 p35.                                          |
| <b>Concentration</b>       | 1 mg/ml                                                                                                               |
| <b>Storage</b>             | Store at -20°C. Avoid repeated freeze/thaw cycles.                                                                    |
| <b>Note</b>                | For research use only                                                                                                 |
| <b>Isotype</b>             | IgG                                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                                            |
| <b>Purity</b>              | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Dilution Range</b>      | WB 1:500 - 1:2000; IHC 1:100 - 1:300; IF 1:200 - 1:1000; ELISA 1:10000;                                               |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                       |