

Product Datasheet

Cleaved-Caspase-9 (D330) antibody [Out of stock] (orb1117440)

Description	Rabbit polyclonal antibody to Cleaved-Caspase-9 (D330)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from the C-terminal region of human Caspase-9 p10.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500 - 1:2000; IHC 1:100 - 1:300; ELISA 1:40000;
Expiration Date	12 months from date of receipt.