

Product Datasheet

OAT3 antibody (orb11178)

Description

OAT3 antibody

Species/Host

Rabbit

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human OAT3. Please contact us for the exact immunogen sequence. The peptide is available as orb374741.

Target

OAT3. OAT3 antibody detects endogenous levels of total OAT3

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

60 kDa

Uniprot ID

[Q9R1U7](#), [O88909](#), [Q8TCC7](#)

NCBI

[001171661](#), [001184732](#)

Entrez


[9376](#)

Dilution Range


IHC-P: 1:50-500, WB: 1:100-2000

Expiration Date


12 months from date of receipt.



IHC-P staining of rat brain tissue using...



Immunohistochemical staining of mouse heart...



IHC-P image of rat kidney tissue using a...