



Product Datasheet

OCT2 antibody (orb11183)

Description	OCT2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from mouse OCT2 (11-100/463 aa)
Target	OCT2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	51 kDa
Uniprot ID	P09086
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:A549
Cell;Primary:Anti-OCT2...



Western blot:Lane
1:Mouse Lymph node;Lan...