


Product Datasheet

OVA antibody (orb11197)

Description	OVA antibody
Species/Host	Rabbit
Reactivity	Gallus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Ovalbumin [hen egg white]
Target	Ovalbumin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43 kDa
CAS Number	4
Uniprot ID	P02768
Dilution Range	WB=1:5000-20000, ELISA=1:5000-100000
Expiration Date	12 months from date of receipt.



Western blot:Ovalbumin protein;Ovalbumin...