



Product Datasheet

P311 antibody (orb11206)

Description	P311 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human P311 (25-68/68 aa)
Target	P311
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	7.9 kDa
Uniprot ID	Q16612
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



IHC-P of human tissue affected by keloid...