



Product Datasheet

MTp53 antibody (orb11211)

Description	MTp53 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr
Immunogen	KLH conjugated synthetic peptide derived from human p53 around the Mutation site N235K N239Y (230-241/393 aa)
Target	Mutant p53 (N235K,N239Y)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21 kDa
Uniprot ID	P04637
Dilution Range	IHC-F=1:100-500, ICC=1:100, Flow-Cyt=1µg/Test, IHC-P=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



IHC-P of human endometrial carcinoma tis...