



Product Datasheet

p63 antibody (orb11214)


Description	p63 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human p63. Please contact us for the exact immunogen sequence. The peptide is available as orb12901.
Target	p63
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	77 kDa
Uniprot ID	O88898 , Q9H3D4 , Q9JJP6
NCBI	001108450 , 001114978
Entrez	8626
Dilution Range	IF/ICC: 1:100-600, IHC-P: 1:100-500, WB: 1:100-1000
Expiration Date	12 months from date of receipt.



IF analysis of rat skin tissue using p63...



IHC-P image of rat stomach tissue using ...



IF image of human lung cancer tissue using p63...