



Product Datasheet

p70 S6 Kinase beta antibody (orb11216)


Description	p70 S6 Kinase beta antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from the middle of human p70 S6 Kinase (265-380/525 aa)
Target	RPS6KB1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	58 kDa
Uniprot ID	P23443
Dilution Range	WB=1:500-2000, Flow-Cyt=1ug/test, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from p...



Immunohistochemical staining of rat hear...



Immunohistochemical staining of rat live...