



Product Datasheet

PACAP 27/38 antibody (orb11220)

Description	PACAP 27/38 antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human PACAP-27/38 (132-169/176 aa)
Target	PACAP-27/38
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	19 kDa
CAS Number	137061-48-4
Uniprot ID	P18509
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



IHC-P of rat brain tissue (dilution at:1...