

Product Datasheet

SARS-CoV-2 Nucleocapsid Protein antibody [Out of stock] (orb1122074)

Description	Rabbit polyclonal antibody to SARS-CoV-2 Nucleocapsid Protein
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	Cell Culture, ELISA, IF, WB
Preservatives	PBS, pH 7.4
Form/Appearance	Liquid
Storage	Store undiluted antibody at 2°C-8°C for up to a week from date of receipt. For long term storage, store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Polyclonal
Purity	Protein A purified
Dilution Range	Optimal dilutions or concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.