

Product Datasheet

SARS-CoV-2 Membrane antibody [Out of stock] (orb1122078)

Description	Rabbit polyclonal antibody to SARS-CoV-2 Membrane
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Preservatives	Contains 0.01% - 0.03 % sodium azide.
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG1
Clonality	Polyclonal
Dilution Range	Optimal dilutions or concentrations should be determined by the end user
Expiration Date	12 months from date of receipt.