

Product Datasheet

SARS-CoV-2 Nucleocapsid antibody [Out of stock] (orb1122084)

Description	Rabbit polyclonal antibody to SARS-CoV-2 Nucleocapsid
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	Cell Culture, ELISA, FC, IF, IHC-P, IP, WB
Preservatives	0.2 µm filtered solution in PBS
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Protein A purified
Dilution Range	Optimal dilutions or concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.