

## Product Datasheet

### SARS-CoV-2 Nucleocapsid antibody [Out of stock] (orb1122085)

|                            |                                                                                                                                                                                                                                                                                |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to SARS-CoV-2 Nucleocapsid                                                                                                                                                                                                                          |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                                                                         |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                                                   |
| <b>Tested Applications</b> | Cell Culture, ELISA, FC, IF, IHC-P, IP, WB                                                                                                                                                                                                                                     |
| <b>Preservatives</b>       | 0.2 µm filtered solution in PBS                                                                                                                                                                                                                                                |
| <b>Form/Appearance</b>     | Liquid                                                                                                                                                                                                                                                                         |
| <b>Storage</b>             | This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt. |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                                          |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                                                                            |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                                     |
| <b>Purity</b>              | Protein A purified                                                                                                                                                                                                                                                             |
| <b>Dilution Range</b>      | Optimal dilutions or concentrations should be determined by the end user.                                                                                                                                                                                                      |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                                                |