

Product Datasheet

SARS-CoV-2 N Protein antibody [Out of stock] (orb1122086)

Description	Rabbit polyclonal antibody to SARS-CoV-2 N Protein
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Preservatives	10 mM PBS, pH 7.2, contains 0.1% sodium azide
Form/Appearance	Liquid
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing.?
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	Optimal dilutions or concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.