

Product Datasheet

SARS-CoV-2 Spike S2 antibody [Out of stock] (orb1122094)

Description	Rabbit polyclonal antibody to SARS-CoV-2 Spike S2
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Preservatives	Citric acid solution
Form/Appearance	Lyophilized, reconstitute with sterile water to a final concentration of 1 mg/ml for use.
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Protein A purified
Dilution Range	Optimal dilutions or concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.