



Product Datasheet

PACAP 38 antibody (orb11221)


Description	PACAP 38 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human PACAP-38 (31-38/38 (162-139/176))
Target	PACAP-38
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	5 kDa
CAS Number	137061-48-4
Uniprot ID	P18509
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, IF=1:100-500
Expiration Date	12 months from date of receipt.



Western blot analysis of rat brain lysat...



rat heart lysates(Lane2) using PACAP 38 ...



IF analysis of rat brain tissue using PA...