

Product Datasheet

HMC Class I (H-2Kb) antibody (orb1122201)

Description	Mouse monoclonal antibody to HMC Class I (H-2Kb)
Species/Host	Mouse
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	Mouse MHC Class I protein H-2Kb
Preservatives	PBS, pH 7.4 with 0.05% sodium azide
Form/Appearance	Liquid
Concentration	100µg/100µl
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Application notes	Optimal dilutions should be calculated by the end user based on their specific application.
Isotype	IgG2a, Kappa
Clonality	Monoclonal
Dilution Range	Optimal dilutions should be calculated by the end user based on their specific application.
Expiration Date	12 months from date of receipt.