

## Product Datasheet

### GFP antibody (orb1122206)

<b>Description</b>	Mouse monoclonal antibody to GFP
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, WB
<b>Immunogen</b>	Full length enhanced green fluorescent protein (aa 1-239)
<b>Preservatives</b>	PBS, pH 7.4 with 0.05% sodium azide
<b>Form/Appearance</b>	Liquid
<b>Concentration</b>	100µg/100µl
<b>Storage</b>	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	Optimal dilutions should be calculated by the end user based on their specific application.
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Dilution Range</b>	Optimal dilutions should be calculated by the end user based on their specific application.
<b>Expiration Date</b>	12 months from date of receipt.