


Product Datasheet


PAF antibody (orb11223)

Description	PAF antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated PAF
Target	PAF
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	0.524 kDa
CAS Number	74389-68-7
Expiration Date	12 months from date of receipt.



IHC-P of rat colon tissue (PAF antibody

...



Immunohistochemical analysis of formalin...