


Product Datasheet

PAK2 (phospho-Ser141) antibody (orb11229)

Description	PAK2 (phospho-Ser141) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PAK2 around the phosphorylation site of Ser141 (YM(p-S)FT)
Target	phospho-PAK2 (Ser141)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23/35/58 kDa
Uniprot ID	Q13177
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse pl...