

Product Datasheet

ROBO1/DUTT1 (Internal) antibody [Out of stock] (orb1123105)

Description	Goat polyclonal antibody to ROBO1/DUTT1 (Internal)
Species/Host	Goat
Reactivity	Bovine, Human, Mouse, Other, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	GDVDLSNKINEMK
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin
Form/Appearance	Liquid
Concentration	100µg/200ul
Storage	This product is stable for several weeks at 4° C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20° C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Affinity purified
Expiration Date	12 months from date of receipt.